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COVID-19 Antigen Rapid Test Clinical Sensitivity and Specificity Study Report

1. Objective

The CLUNGENE[®] COVID-19 Antigen Rapid Test (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

2. Method

A study of 770 direct nasopharyngeal swabs was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasopharyngeal swabs were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasopharyngeal swab was tested directly with the COVID-19 Antigen Rapid Test according to product instructions for use, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Swabs were randomly assigned to testing with the COVID-19 Antigen Rapid Test or RT-PCR assay and were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$. The negative percent agreement (NPA) was calculated as $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$. Accuracy was calculated as $100\% \times ([\text{True Positive} + \text{True Negative}] / \text{Total sample Qty})$. The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene β -Actin is not higher than 35.

4. Enrollment criteria (inclusion/exclusion criteria)

4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The

sensitivity was calculated for the different Ct value range (Ct value \leq 33 and Ct value \leq 37).

COVID-19Antigen		RT-PCR (Ct value \leq 33)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct \leq 33): 98.0% (145/148), (95%CI: 94.2% ~99.3%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~99.9%)

Accuracy: 99.3% ((145+593)/743), (95%CI: 98.4% ~99.7%)

COVID-19Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct \leq 37): 92.0% (161/175), (95%CI: 87.0% ~95.2%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~99.9%)

Accuracy: 97.9% ((161+593)/770), (95%CI: 96.6% ~98.7%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

6. Conclusion

Taken together, the CLUNGENE[®] COVID-19 Antigen Rapid Test had a positive percent agreement (sensitivity) of 98.0% (95%CI: 94.2% ~99.3%) with specimens of a Ct count \leq 33, 92.0% (95%CI: 87.0% ~95.2%) with specimens of a Ct count \leq 37, negative percent agreement (specificity) of 99.7% (95%CI: 98.8% ~99.9%), and accuracy of 99.3% (95%CI: 98.4% ~99.7%) with specimens of a Ct count \leq 33, 97.9% (95%CI: 96.6% ~98.7%) with specimens of a Ct count \leq 37.

COVID-19 Antigen Rapid Test Cassette (Nasal Swab) Clinical Sensitivity and Specificity Study Report

1. Objective

The CLUNGENE[®] COVID-19 Antigen Rapid Test Cassette (Nasal Swab) (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

2. Method

A study of 617 direct nasal swabs was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasal swabs were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasal swab was tested directly with the COVID-19 Antigen Rapid Test Cassette (Nasal Swab) according to product instructions for use, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Swabs were randomly assigned to testing with the COVID-19 Antigen Rapid Test Cassette (Nasal Swab) or RT-PCR assay and were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$. The negative percent agreement (NPA) was calculated as $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$. Accuracy was calculated as $100\% \times ([\text{True Positive} + \text{True Negative}] / \text{Total sample Qty})$. The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene β -Actin is not higher than 35.

4. Enrollment criteria (inclusion/exclusion criteria)

4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value \leq 30 and Ct value \leq 37).

COVID-19 Antigen		RT-PCR (Ct value \leq 30)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	117	3	120
	Negative	3	462	465
Total		120	465	585

PPA (Ct \leq 30):97.5% (117/120), (95%CI: 92.9%~99.2%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

COVID-19 Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct \leq 37):91.4% (139/152), (95%CI: 85.9%~94.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

6. Conclusion

Taken together, the CLUNGENE[®] COVID-19 Antigen Rapid Test Cassette (Nasal Swab) had a positive percent agreement (sensitivity) of 97.5% (95%CI: 92.9%~99.2%) with specimens of a Ct count \leq 30, 91.4% (95%CI: 85.9%~94.9%) with specimens of a Ct count \leq 37, negative percent agreement (specificity) of 99.4% (95%CI: 98.1%~99.8%).

EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Clongene Biotech Co., Ltd.
No.1 Yichuang Road, Yuhang Sub-district
Yuhang District
311121 Hangzhou
China**

We declare under our sole responsibility that

the medical device: **COVID-19 Antigen Rapid Test**

of class: **Other**
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III**

Applicable standards: **EN ISO 13485:2016 EN ISO 15223-1:2016
EN ISO 23640:2015 EN13612:2002/AC:2002
EN 13975:2003 EN ISO 14971:2012
EN ISO 18113-1:2011 EN ISO 18113-2:2011
EN 62366-1:2015**

Name and address of the authorized representative: **Shanghai International Holding Corporation GmbH (Europe)
Eiffestrasse 80
20537 Hamburg
Germany**

Hangzhou, July.15.2020

Place, date

Shujian Zheng, Legal representative

Name and function

COVID-19 Antigen Rapid Test

The COVID-19 Antigen Rapid Test is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab, nasal swab or oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

Features

CE

Easy to collect samples.

No equipment required.

Instant result at 15 minutes.

Results are clearly visible.


Suitable for large-scale rapid screening.

Components

Extraction Reagent
Tube x 25

Sterilized Swab
x 25

Test Cassette
x 25

Work Station
x 1


Package Insert
x 1


Samples

Nasopharyngeal Swab


Nasal Swab

Oropharyngeal Swab

Test Procedure

Performance Characteristics

Clinical Performance

Nasopharyngeal Swab

770 nasopharyngeal swabs were collected from individual symptomatic patients by using CLUNGENE® COVID-19 Antigen Rapid Test. The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value ≤ 33 and Ct value ≤ 37).

Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value ≤ 33)		Total
		Positive	Negative	
CLUNGENE®	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct ≤ 33): 98.0% (145/148), (95%CI: 94.2%~99.3%)
 NPA: 99.7% (593/595), (95%CI: 98.8%~99.9%)

COVID-19 Antigen		RT-PCR (Ct value ≤ 37)		Total
		Positive	Negative	
CLUNGENE®	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct ≤ 37): 92.0% (161/175), (95%CI: 87.0%~95.2%)
 NPA: 99.7% (593/595), (95%CI: 98.8%~99.9%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

Nasal Swab

617 nasal swabs were collected from individual symptomatic patients by using CLUNGENE® COVID-19 Antigen Rapid Test. The sensitivity was calculated for the different Ct value range (Ct value \leq 33 and Ct value \leq 37).

Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value \leq 33)		Total
		Positive	Negative	
CLUNGENE®	Positive	132	3	135
	Negative	4	462	466
Total		136	465	601

PPA (Ct \leq 33): 97.1% (132/136), (95%CI: 92.7%~98.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

COVID-19 Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
CLUNGENE®	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct \leq 37): 91.4% (139/152), (95%CI: 85.9%~94.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

Limit of Detection (Analytical Sensitivity)

The study used cultured SARS-CoV-2 virus which is heat inactivated and spiked into nasopharyngeal swab specimen. The Limit of Detection (LoD) is 5.7×10^2 TCID₅₀/mL.

Cross Reactivity (Analytical Specificity)

32 commensal and pathogenic microorganisms that may be present in the nasal cavity have been evaluated and no cross-reactivity was observed.

Interference

17 potential interference substances with different concentration were evaluated and found no affect to the test performance.

High-dose Hook Effect

The COVID-19 Antigen Rapid Test was tested up to $1.0 \times 10^{5.67}$ TCID₅₀/mL of inactivated SARS-CoV-2 and no high-dose hook effect was observed.

Hangzhou Clongene Biotech Co., Ltd.

Add: No.1 Yichuang Road, Yuhang Sub-district, Yuhang District, 311121 Hangzhou, China

Tel: +86 (0)571-88262120

Fax: +86 (0)571-88261752

EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Clongene Biotech Co., Ltd.**
No.1 Yichuang Road, Yuhang Sub-district
Yuhang District
311121 Hangzhou
China

We declare under our sole responsibility that

the medical device: **COVID-19 Antigen Rapid Test Cassette (Nasal Swab)**

of class: **Other**
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III, excluding Section 6**

Applicable standards: **EN ISO 13485:2016** **EN ISO 15223-1:2016**
EN ISO 18113-1:2011 **EN ISO 18113-2:2011**
EN ISO 23640:2015 **EN ISO 14971:2019**
EN 13641:2002 **EN13612:2002/AC:2002**
EN 13975:2003 **EN 62366-1:2015**
ISO 17511:2020

Name and address of the authorized representative: **Shanghai International Holding Corporation GmbH (Europe)**
Eiffestrasse 80
20537 Hamburg
Germany

Hangzhou, February 04, 2021.

Place, date

Shujian Zheng, Legal representative

Name and function

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
**Hangzhou Clongene Biotech
Co., Ltd.**
**No. 1 Yichuang Road, Yuhang Sub-district
Yuhang District
311121 Hangzhou
P.R. China**

has established and applies a quality management system for medical devices
for the following scope:

**Design and Development, Manufacture and Distribution of
In-vitro Diagnostic Rapid Test of Fertility, Drug of Abuse,
Infectious Diseases, Tumour Markers and Cardiac Markers**

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-11-25
Certificate Registration No.: SX 60152722 0001
An audit was performed. Report No.: 15073650 008
This Certificate is valid until: 2023-11-18

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July 6, 2021

To whom it may concern,

We, as the manufacturer of the CLUNGENE COVID-19 Antigen Rapid Test, assure that the following major SARS-CoV-2 variants which seem to spread more easily and quickly than other variants and are circulating globally, can be detected by CLUNGENE COVID-19 Antigen Rapid Test:

- ✓ **B.1.1.7 (Alpha):** This variant was initially detected in the UK;
- ✓ **B.1.351 (Beta):** This variant was initially detected in South Africa in December 2020;
- ✓ **P.1 (Gamma):** This variant was initially identified in travelers from Brazil, who were tested during routine screening at an airport in Japan, in early January;
- ✓ **B.1.617.2 (Delta):** This variant was initially detected in India.

The recombinant N proteins of these variants can be detected by the CLUNGENE COVID-19 Antigen Rapid Test, at concentrations same as that of the ordinary recombinant N protein without variation (amino acid: Met1-Ala419, accession # QHO62884.1).

Yours sincerely,

Shujian Zheng,

Managing Director of Hangzhou Clongene Biotech Co., Ltd.

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EC REP Shanghai International Holding Corp.GmbH (Europe)
Pitea S.r.l.s