

## Declaration on the detection validity of omicron variant

To Who It May Concern,

For the newly discovered variant of novel coronavirus Omicron (B.1.1.529), we have analyzed WIZ SARS-CoV-2 Antigen Rapid Test Kit and came to a conclusion: there is no off-target or missing detection occurred in the latest test for variant of Omicron (B.1.1.529). WIZ SARS-CoV-2 antigen rapid test adopts double antibody sandwich method to detect Nucleocapsid protein of SARS-CoV-2 viruses, Omicron NP protein mutation sites were P13L, 31-33, R203K, G204R. WIZ SARS-CoV-2 antigen rapid test detectable sites are not at these mutation sites. Our team will continue to track the latest variant of novel coronavirus to ensure stable supply.



<u>Tel:+86-592-6808278</u> Fax: +86-592-6808279

www.wizbiotech.com